



Options & Accessories

Eurus Options & Accessories

Chair



Double articulating headrest



Power headrest




Upholstery Luxury type

SH vinyl colour list for luxury upholstery


When you select the luxury upholstery option, the colour of vinyl should be chosen from the following SH vinyl colour list.



SH9
Dark Blue



SH1
Black



SH7
Maroon



SH18
England Green



SH22
Orange



SH10
Camel



SH6
Mustard Yellow



SH15
Ivory

Standard ● Selectable △ Not available —

Chair	S1 & S3 Contour	S4 Folding legrest
Double articulating headrest (2 integrated stick switch, manual & auto mode chair switch)	●	△
Power headrest (3 integrated stick switch, manual, auto mode chair switch & power headrest)	△	●
Armrest at left side only (fixed)	●	—
Armrest at both side, Left side fixed & right side rotated	△	●
Upholstery seamless type	●	●
Upholstery seamless soft type	△	—
Upholstery luxury type	△	△

Dentist Element




Place type



Holder type



Belmont Rod system



Luzzani 6way-Syringe



USB-Port



L-CAM HD SHOT (Straight type)

Standard ● Selectable △ Not available —


Dentist Element	Holder type	Belmont Rod system
Syringe		
Luzzani 3-way syringe	△	△
Luzzani 6-way syringe	△	△
Air turbine hose with light		
Midwest 4-way turbine hose with fiber optic	△	△
Micro motor (*1)		
Bien Air MCX 1pc	△	△
Bien Air MCX 2pcs	△	△
Bien Air MCX 3pcs	△	△
Bien Air MX2 1pc	△	△
Bien Air MX2 2pcs	△	△
Bien Air MX2 3pcs	△	△
NSK NLX plus 1pc	△	△
NSK NLX plus 2pcs	△	△
Ultrasonic scaler		
SATELEC SP 4055 Newtron	△	△
SATELEC Xinetic	— (*2)	— (*2)
NSK Varios 170	△	△
EMS Piezon No Pain	— (*2)	— (*2)
Intraoral-Camera		
L-CAM HD SHOT (Straight type or angled type)	△	△
ACTEON Sopro 617	△	△
ACTEON Sopro Care	△	△
HP hose with quick connector for dentist and assistant element (5+1)	△	△
Small size tray	△	△
Large size tray	△	△
Water heater for handpieces	△	△
Intraoral film viewer	△	△
Panorama film viewer	△	—
Removable holder sleeve	●	—
Removable handpiece pads	—	●
USB Port	△	△
Place type handpiece holder	△	—

*1: When you order 2 or 3 micro motors, you can't select the two from different kinds of micro motors. (i.e. one Bien Air MCX and one Bien Air MX2) In this case, you need to select them from one kind.

*2: Please contact us for the availability.

Eurus Options & Accessories

Unit Body



Vacuum cleaning system (Dürr Hygosuc)



Sensored cupfiller & bowl Flush



Clean water system with 2L water tank



Separated junction box



LED indicator



Wireless foot controller



Monitor Bracket (not incl. monitor)



Built-in flush-out adapters

Standard ● Selectable △ Not available -

Unit body	Chair mount	Pedestal
Electric timered cupfiller & bowl flush	●	●
Sensored cupfiller & bowl flush	△	△
Rotating cuspidor bowl (porcelain)	●	●
Water heater for cupfiller	△	△
Clean water system		
Clean water system with 2L water tank (FF-version, replacing air gap)	△	△
Clean water system with 2L water tank (SF-version, not replacing air gap)	△	△
Metasys WEK (Water decontamination system with air gap)	△	△
Metasys WEK Light (Water decontamination system without air gap)	△	△
Separator		
Dürr Spittoon bowl valve	△	△
Dürr CS1 (Separator) & Rinsing system (*3)	△	△
Dürr CAS1 (Amalgam collector) & Rinsing system (*3)	△	△
Metasys Eco1 (Separator)	△	△
Metasys MST1 (Amalgam collector)	△	△
Built-in flush-out adapters	- (*4)	- (*4)
Vacuum cleaning system (Dürr Hygosuc)	-	- (*4)
LED indicator	●	●
Wireless foot controller	△	△
Extra switch	△	△
Monitor bracket	△	△
Neovo Monitor 22 inch white with handle (MPG-conformity with glass protection plate, connectors: DVI/VGA/HDMI/S-VHS)	△	△
Service outlet for Water	△	△
Service outlet for Air	△	△
Separated Junction Box	△	△

*3 : You can't select the combination of Dürr CS1, CAS1 and Metasys WEK.

*4 : Please contact us for the availability.

Assistant Element



Ratchet gear type
assistant arm

Standard ● Selectable △ Not available -

Assistant-Element	Chair mount	Pedestal
Ratchet gear type assistant arm	△	△
Luzzani 3-way syringe	△	△
HP hose with quick connector for assistant syringe (see dentist element)	△	△
Dürr high volume ejector with ball joint & Saliva ejector handpiece (Dürr)	●	●
Twin vacuum system	●	●
Vacuum hose clamp	△	△

Eurus Options & Accessories

Light



Patient mirror

Stool



Eurus Stool

Standard ● Selectable △ Not available -

Light

Patient mirror	△
----------------	---

Stool

Eurus stool	△
-------------	---

Dimensions









S1



S3



S4



MH:Double articulating headrest
PH: Power headrest

